


## Use of Microfluidizer™ technology for graphene exfoliation

This Application Note gives an overview of the properties and applications of graphene. It also discusses preparation techniques and a case study executed in collaboration with the University of Cambridge to produce graphene for use as conductive inks. [1]

### What is graphene?

Graphene is a material comprised of a sheet of carbon just one atom thick. See structure below:



Graphene has many useful properties:

- It is ultra-light yet extremely tough.
- It is 200 times stronger than steel, but it is incredibly flexible.
- It is the thinnest material possible as well as being transparent.
- It is a superb conductor.
- It can act as a perfect barrier - not even helium can pass through it.


These properties lead to advances in many potential applications such as:

- Membranes to clean drinking water
- High strength composites and coatings
- High quality batteries.
- Microfluidic-based sensors
- Many electronic applications

There have been many different methods investigated for production of graphene. The first method used was the use of adhesive tape to exfoliate graphene from graphite. There have been many other techniques discovered including chemical vapor deposition (CVD),  $\text{CO}_2$  reduction, spin coating and other complex processes. These processes have all been developed on the research scale and would be difficult/expensive to use for bulk production.

- Batch variation or probe variation can only spread concentration of up to ~0.2 g/L, with yield of only 1-2%.
- High speed rotor-stator mixer gave much more results of concentration <math>< 1 \text{ g/L}</math>, and yield <math>< 10\%</math>.
- With Microfluidizer, a concentration of ~100 g/L with 100% yield was achieved. This is 500-1000 times higher in terms of product concentration and 50-100 times higher in terms of yield.

The Microfluidizer process can generate shear forces of  $10^7 \text{ s}^{-2}$  which liberates the graphene sheets from primary graphite (see image below). This process is encouraging because Microfluidizer technology is already scalable to hundreds of liters per hour meaning that it can be used for large scale manufacturing of graphene.



SEM images taken from castings comparing (a) starting graphite, (b) after 1 cycle, (c) after 5 cycles, and (d) after 200 cycles through the Microfluidizer processor at 20,000 psi. The scale bar is 5  $\mu\text{m}$ .