

## PMX 220 ZetaView TWIN Laser

Standard Technical Data (availability depending on selected modules)

### General Features

Measurement Principle:	<ul style="list-style-type: none"> <li>● Precision-engineered motorized scanning Nanoparticle Tracking Analysis (NTA) instrument for tracking the movement of individual visualized nanoparticles in suspension</li> <li>● Real-time visualization of Brownian Motion and electrophoretic mobility, for measuring size, concentration and zeta potential in scattering and fluorescence modes.</li> <li>● Two simultaneous aligned and software-controlled lasers for fluorescence measurements</li> <li>● Software controlled emission filter wheel for quick changes between the fluorescence measurements</li> <li>● Fast scanning to acquire and analyze typically 1000 particles in ~ 1 minute</li> <li>● Software-controlled pumps for flushing and sample sub-dosing</li> </ul>
Samples:	<ul style="list-style-type: none"> <li>● Nanoparticles suspended in polar liquids (e.g. water, alcohols) for size, concentration, fluorescence and zeta potential studies</li> </ul>

### Hardware

Equipment:	<ul style="list-style-type: none"> <li>● ZetaView® PMX-220 TWIN Laser main unit is equipped with Cell Assembly, two simultaneous aligned lasers (see section Lasers) and bottles for buffer rinse</li> <li>● Power of statistics by automated unique scan and dose control for measurement of 1 - 100 independent sub-volumes</li> <li>● Zeta potential option</li> <li>● Software controlled double Fluorescence option features short acquisition times to avoid negative effect of photo bleaching</li> </ul>
Optical Layout :	<ul style="list-style-type: none"> <li>● 90° laser scattering video microscope with x10 magnification</li> <li>● Automated focusing of laser and microscope</li> </ul>
Laser sets:	<ul style="list-style-type: none"> <li>● Special TWIN Laser design</li> <li>● Available laser wavelengths combinations: <ul style="list-style-type: none"> <li>405 nm / 488 nm</li> <li>405 nm / 520 nm (coming soon)</li> <li>405 nm / 660 nm (coming soon)</li> <li>488 nm / 660 nm (Available approx. Q3/2018)</li> </ul> </li> <li>at typical laser power of &gt;30 mW per laser</li> <li>● Pulse duration each laser 0.1 ms up to continuous</li> </ul>
Camera:	<ul style="list-style-type: none"> <li>● Sensitive CMOS camera 640 x 480 pixels</li> <li>● Variable frame rate from 1 to 60 Hz for optimum resolution and fast acquisition</li> </ul>

# PMX 220 ZetaView TWIN Laser

## Standard Technical Data (availability depending on selected modules)

<b>Fluorescence Filters:</b>	<ul style="list-style-type: none"> <li>● Software controlled, automated filter wheel</li> <li>● Available long wave-pass (LWP) filter combinations: <ul style="list-style-type: none"> <li>430 nm / 500 nm for 405 / 488 laser combination</li> <li>430 nm / 550 nm for 405 / 520 laser combination (coming soon)</li> <li>430 nm / 680 nm for 405 / 660 laser combination (coming soon)</li> <li>500 nm / 660 nm for 488 / 660 laser combination (Available approx. Q3/2018)</li> </ul> </li> <li>● Bandpass filter available on request</li> </ul>
<b>Cell Assembly:</b>	<ul style="list-style-type: none"> <li>● Z-NTA – slide-in assembly for size, concentration and dual fluorescence measurements plus zeta potential experiments in aqueous and organic solvents with pumps for 2 different liquids/buffers – for rinsing and sub-dosing experiments, electrical field sensing</li> </ul>
<b>Cleaning:</b>	<ul style="list-style-type: none"> <li>● Cell cleaning recommended weekly – cell resistant to &gt;1000 brush cleanings</li> <li>● Cleaning of driver electrodes required after more than 1000 zeta potential runs</li> <li>● Cleaning kit and basic replacement parts included in delivery</li> </ul>
<b>Temperature Range/Control:</b>	<ul style="list-style-type: none"> <li>● Working external temperature range: 5°C to 45°C</li> <li>● Sample temperature control: Peltier temperature control from RTP-5oC to 55oC with dew-point sensing</li> </ul>

### Software

<b>Communication:</b>	<ul style="list-style-type: none"> <li>● Software provided on pre-configured PC, communication via Ethernet</li> </ul>
<b>Quality Control:</b>	<ul style="list-style-type: none"> <li>● Cell quality check, daily performance check, outlier control with automatic Grubbs statistical analysis of measurement data</li> </ul>
<b>Live Monitoring:</b>	<ul style="list-style-type: none"> <li>● Number of detected particles in scatter or fluorescence, scattering intensity, conductivity*, temperature, particle drift</li> </ul>
<b>Standard Operating Procedures (SOP):</b>	<ul style="list-style-type: none"> <li>● Fully customisable SOPs for different samples/applications</li> </ul>
<b>Analysis and Reports:</b>	<ul style="list-style-type: none"> <li>● Data Analysis in scatter and / or fluorescence mode: particle size distribution profiles, concentration, overlays and averaging, scatter plots, zeta-potential distribution profiles, sub-population analysis</li> <li>● Data export format: AVI, TXT, CSV, FCS</li> <li>● PDF reports containing key results</li> </ul>

### Measurement Specifications

<b>Size/ Concentration:</b>	<ul style="list-style-type: none"> <li>● Concentration range: 10<sup>5</sup> – 10<sup>9</sup> particles/ml</li> <li>● Particle size: 15nm – 1000nm (dependent on sample and laser selection)</li> <li>● Accuracy: ±5nm (for 100nm polystyrene latex)</li> <li>● Reproducibility: ±2nm (for 100nm polystyrene latex)</li> </ul>
<b>Fluorescence:</b>	<ul style="list-style-type: none"> <li>● Concentration range: 10<sup>5</sup> – 10<sup>9</sup> particles/ml</li> <li>● Particle size: 30nm – 1000nm (dependent on fluorescent dye and laser selection)</li> <li>● Accuracy: ±5nm (for 100nm polystyrene latex)</li> <li>● Reproducibility: ±2nm (for 100nm polystyrene latex)</li> </ul>

# PMX 220 ZetaView TWIN Laser

## Standard Technical Data (availability depending on selected modules)

Zeta Potential*:	<ul style="list-style-type: none"><li>● Working range: -500 to +500mV</li><li>● Concentration range: <math>10^6</math> – <math>10^{10}</math> particles/ml</li><li>● Particle size: 20nm – 5000nm (dependent on sample and laser selection)</li><li>● Conductivity range: 3<math>\mu</math>S/cm – 15mS/cm</li><li>● Accuracy: <math>\pm 4</math>mV (for alumina zeta potential standard)</li><li>● Reproducibility: <math>\pm 2</math>mV (for alumina zeta-potential standard)</li></ul>
General:	<ul style="list-style-type: none"><li>● Minimum sample quantity: 500<math>\mu</math>l of sample at <math>10^5</math> particles/ml</li><li>● pH range: 2 – 12</li><li>● Temperature: 5°C to 45°C (external temperature)</li><li>● Sample volume visualised and tracked by the camera for a single measurement: 11 x 3nL</li></ul>
Reference Materials:	<ul style="list-style-type: none"><li>● Nominal 100 nm reference suspension for size</li><li>● Two nominal 100 nm reference suspensions for fluorescence</li><li>● Nominal +50mV reference suspension for zeta potential*</li></ul>

### Dimensions.

Physical:	<ul style="list-style-type: none"><li>● Footprint (W x D x H): 20 x 30 x 25cm</li><li>● Weight: 8.5kg (main unit, PC extra)</li><li>● Shipping box with standard content: 48 x 62 x 63cm; 22kg</li></ul>
Electrical:	<ul style="list-style-type: none"><li>● 90-240V, 47-63Hz, 50VA</li></ul>

### Warranty & Support

Warranty:	<ul style="list-style-type: none"><li>● 1 year (glass excluded).</li></ul>
Service & Support:	<ul style="list-style-type: none"><li>● Reaction time: 48 hours</li><li>● Maintenance, service and IQ/OQ/PQ contracts available on demand#</li><li>● Support via telephone, email and TeamViewer for trained users free of charge during warranty period</li><li>● Training courses for new users available on demand</li><li>● Special arrangements and specifications available on demand – quotation required</li></ul>

\* With 'Z-NTA' cell assembly only.

Version ZV220 EN 1.0 April 2018



Manufactured by Particle Metrix, distributed in the UK and Ireland by **analytik**.