



HYPERAI|XPERT

Multi-Sensor Laboratory
Phenotyping System

LemnaTec Alxpert Product Family

Documentation, Phenotyping, Seed Testing and more

HYPERAI|XPERT

The most flexible multi-sensor laboratory phenotyping system for Arabidopsis, seedlings, petri dishes, MTPs, and many more sample types.

- Rapid digital analysis of phenotypic traits
- Large set of sensor options
- Advanced analytical software and machine learning
- Standardized and repeatable data acquisition
- High throughput
- Pots, trays, dishes, MTPs, beakers
- Adaptors for many sample types
- Up to 40 cm sample height with manual loading.
Up to 20 cm with automated sample loader (TrayProvider™)

Features

- Experiment-driven design
- High-resolution phenotyping
- Comprehensive machine-learning based analytical software
- All major phenotyping cameras and sensors available
- Configurable multi-sensor imaging
- Automatic sample loading available
- Integrable into climate chambers

Multisensor imaging

- Visible light images provide information on sample dimensions, morphology, and color
- An IMAGING-PAM camera delivers dynamic chlorophyll fluorescence parameters.
- Hyperspectral imaging delivers a range of wavelength data for calculation of known and user-defined vegetative indices.
- Near-infrared (NIR) images capture the intensity of reflected NIR radiation at 1450 nm to measure water content.
- Fluorescence images capture the intensity of a wide range of fluorescent pigments.
- 3D laser scanning generates point clouds for detailed architectural and morphological information.

Configuration and options

- The Base Module supports manual input of microtiter plates, petri dishes, plants in pots or plant flats.
- Cameras and sensors move across the sample stage inside the base module and are positioned with high accuracy
- The Sensor Exchange System stores cameras and sensors and the moving arm attaches them for running image acquisitions across the sample stage
- Tray Adaptors for use with trays, MTPs and other sample carriers
- Light table for MTPs, petri dishes etc. available
- The TrayProvider™, an automatic tray loading module enables high throughput sample handling for efficient pipeline screening.
- The HyperAlxpert can be integrated into Climate Chambers



individual plant pots
and plant trays



Sample trays
(here: seeds)



Microtiter plates
(MTPs)

HyperAlxpert provides various possibilities in research and development.

GENETICS

screening of populations, assessment of gene functions,
QTL discovery, biomarker studies

PLANT CULTIVATION PRODUCTS

compound screening, treatment efficiency

PLANT RESEARCH AND BREEDING

genotype-phenotype links, environmental effects, physio-
logical effects


PLANT STRESS AND DISEASES

resistance responses, disease studies


Application examples



Population screen with seedlings
– genotype or treatment effects
Growth studies
Phenotypic assessments
Stress response monitoring



Chlorophyll Fluorescence studies
Status and activity of photosystem
Dynamic chlorophyll fluorescence imaging




Plant disease rating and plant health assessments
Physiological phenotyping




Fluorescent biomarkers in gene expression studies


Assessment of samples in petri dishes




Compound or genotype screen in MTPs



Plant performance, development, and stress studies



3D imaging of plant architecture



Manufactured by LemnaTec, distributed in the UK and Ireland by **analytik**.